



Checklist For Pre-analytical Sample Handling

Patient's code: [redacted] Center's code: [redacted] Investigator's code : [redacted]
Matrix-box code: [redacted]

Patient has given written informed consent.	
Patient in the fasting state.	
Spinal tap performed in a sitting position.	
First approach of spinal tap successful.	
CSF transparent on visual inspection.	
Blood, serum and plasma collected the same day as CSF.	
The spinal tap was finished at o'clock. The CSF was aliquoted at o'clock.	
CSF preparation completed within 30 min.	
Blood, serum, and plasma processed within 1.5 h after collection.	
Positions in the Matrix-box : 1-6.	
Positions in the Matrix-box : 7-12.	
The material has never been thawed after being frozen.	
The complete routine CSF and serum work out is attached.	

Comments: